

# i-Glide *and* i-Ject

## *Cartridges and Injectors*

### **i-Glide disposable Cartridges & i-Ject disposable Injectors.**

i-Medical® "i-Glide" cartridge, used in conjunction with our "i-Ject" disposable injector, has a hydraulic effect on the viscoelastic solution, driving the hydrophilic one piece IOL with a silicone stamp. Thinner and smaller than the traditional cartridge, it ensures the safety of the patient's eye and the trueness of the Intraocular Lens during injection.

Because it is thinner and smaller, the surgeon, through the IOL injection technique, can make a micro incision of only 1.85 mm - 2.5 mm. A wonderful improvement on previous methods.




**Typ / Order-No:** i-Glide 1.85 Small

as a Set incl. i-Ject Injector and Silicone Stamp

#### **Theoretical technical information:**

- Oval Shape
- Outer: 1.85 mm x 2.1 mm
- Incision > 1.85 / <2.5 mm



i-Ject Injector including Silicone Stamp.